

2012		EPOND accessories for GV Instruments IsoProbe and Platform ICP-MS
EPOND P/N°	OEM P/No.	Part Descriptions
2050-0700N		Standard Concentric Glass Nebulizer, with 'NebLink'
2050-0700NQ		Standard Concentric Quartz Nebulizer, with 'NebLink'
2050-0710NT		High-Performance 'Tyfoon' Concentric Glass Nebulizer
2050-M030-C1	TR30-C1	<b>'Meinhard'</b> Concentric glass nebulizer, Type C
2050-4252		<b>'Boreal'</b> 'NB-Class' ( No-Break ) Concentric Polymer Nebulizer
2550-2020		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 20 ul/min uptake version
2550-2000		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 50 ul/min uptake version
2550-2001		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 75 ul/min uptake version
2550-2002		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 100 ul/min uptake version
2550-2003		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 200 ul/min uptake version
2550-2005		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		<b>PolyPro-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		<b>PFA-ST</b> HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2045		Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min.
2550-2043		Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min.
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
		<b>Capillaries with integrated probes - available for all common autosamplers</b>
2045-3305		Standard <b>'Electron'</b> 22ml. Mini-cyclonic glass spray chamber
2045-3307		Standard <b>'Positron'</b> baffled Mini-cyclonic glass spray chamber
6609-5853	6368315	Jacketed <b>'Electron'</b> 22ml. Mini-cyclonic glass spray chamber
2045-3305Q		Quartz version of 'Electron' Mini-cyclonic spray chamber
2045-3301		Standard <b>'Althea'</b> 50ml. Cyclonic glass spray chamber
2045-3303		Standard <b>'Libra'</b> 50ml. Baffled Cyclonic glass spray chamber
2045-3301Q		Quartz version of 'Althea' 50ml. Cyclonic spray chamber
2045-3321	6368314	Jacketed <b>'Althea'</b> 50ml. Cyclonic glass spray chamber
EP-4157-1000	ES-4157-1000	<b>PC3</b> Peltier-cooled ESI Quartz Spray Chamber, with Power Supply
EP-9100-5400	ES-9100-5400	<b>PC3-LT</b> (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2045-0220		VG Jacketed Spray Chamber with Impact Sphere
2045-022Q		VG Quartz Jacketed Spray Chamber, with Impact Sphere
2045-0130		VG Quartz Conical Spray Chamber, for Peltier system
2045-0150		Glass Coaxial Spray Chamber, Water Cooled
2045-0160		Quartz Coaxial VG Spray Chamber, Water Cooled
2045-0800		HF-resistant Teflon Baffled Spray Chamber
2045-0804		HF-resistant Ultem Cyclonic Spray Chamber
2045-0050		Standard Quartz Torch 1.5mm Injector
2045-0060		Standard Quartz Torch 2mm Injector
2045-0060		Extended Quartz Torch 1.5mm Injector
2045-0080		Demountable Torch Body, for 2045-0100
2045-0100		Semi-Demountable Torch, with Qtz Injector
EP-1202-1802		Semi-Demountable Torch, O-ring free, with Sapphire Injector and spray chamber adaptor
EP-1203-0182		Sapphire injector for O-ring free Semi-Demountable Torch, with spray chamber adaptor
EP-1202-0002		O-ring free Semi-Demountable Torch with PFA base, without injector
2045-0090		Quartz Injector, for 2045-0100
2045-0110		Alumina Injector, for 2045-0100
2045-0250		Standard VG Torch Bonnet, for PQ series, Micromass
2045-0310		VG Elbow, Ball/Socket Pyrex
2045-0320		VG Elbow, Ball/Socket Quartz

EPOND No.	OEM P/No.	Micromass Isoprobe ( GV Instruments )
M1001-Ni/Cu	M695080BD1	Ni Sampler Cone with Cu base
M1002-Ni	M695011BD1	Nickel Skimmer Cone
M1006-Pt		Platinum Sampler Cone, 15.2mm insert
M1008-Pt		Platinum Skimmer Cone, 7.6mm insert
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe
EPOND No.	OEM P/No.	New Isoprobe and Platform cones
MM1001-Ni	M749629CD3	Nickel sampler cone
MM1001-AI	M749629CD2	Aluminium Sampler Cone
MM1002-Ni	M695440BD1	Nickel Skimmer Cone
MM1002-AI	M695440BD2	Aluminium Skimmer Cone
MM1002X-Ni		Nickel Skimmer with enhanced low-mass performance
MM1002X-AI		Aluminium Skimmer with enhanced low-mass performance
MM3000		Copper cartridge
MM5003	M749634BC3	QuickCone' set of cones - ( MM1001-Ni + MM1002-Ni + MM3000 cartridge )
MM1003-Ni	M749666B	Ni collimator cone
MM1008-Pt		Platinum Skimmer Cone
		<i>MM' parts require an MM3000 copper cartridge for mounting, suit both types ICP-MS</i>
EPOND No.	OEM P/No.	Micromass Platform
MP1001-Ni/Cu	M695439BD1	Ni Sampler Cone with Cu base
MP1001-AI	M695439BD2	Aluminium Sampler Cone
MP1002-Ni	M695440BD1	Nickel Skimmer Cone
MP1002-AI	M695440BD2	Aluminium Skimmer Cone
MP1006-Pt/Cu		Platinum Sampler Cone, copper base
MP1008-Pt		Platinum Skimmer Cone
EPOND No.		Accessory Description
APEX-Q-VG		<b>APEX-Q</b> Desolvator, with 2 PFA concentric nebs and 'Pluto' flow monitor
2019-4260		<b>'XENA'</b> micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
8100-2051		<b>'PEARL'</b> Cleaning Device for concentric glass nebulizers
ACM-E2		ESI Nafion <b>Active Membrane Device</b> , with DC Power Supply
EP-4599-3002		<b>Spiro TMD</b> Membrane Desolvation Module ( latest model of ACM ), with Power Supply
6808-2280		<b>'ISOBAR'</b> Internal Standard Kit, HF-resistant, with tubing and connectors
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		<b>EPOND 'RubeTubes'</b> , sample tube connectors, packet of 40
DLK-075		<b>DrainLink</b> nebulizer adaptor and tube with 0,75mm ID, pkt 5
NLK-030		<b>Neb-Link</b> sample tubes, packet of 10, ID=0,30 mm, L=60cm.
<b>Autosamplers and ICP / AAS Standard and Reference Solutions</b>		
SC-8000-XXXX		Complete range of ESI SC-series autosamplers - please ask for details
ICP-standards		Certified ICP / AAS Standard solutions - <a href="http://www.e-pond.biz/standards">http://www.e-pond.biz/standards</a>
		<b>we do custom standards too ! ask for a quote ...</b>
<b><a href="http://www.e-pond.biz">www.e-pond.biz</a></b>		