

2012		EPOND Parts List, PE Optima 2000DV – 8000DV
EPOND N°	OEM No.	Descriptions
2050-0700N	N810-2011	Standard ' Caliber ' Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-0700NQ		Standard Concentric Quartz Nebulizer, with ' <i>NebLink</i> ' sample tube
2050-0710N	N068-1574	Standard Concentric Glass Nebulizer, with <i>NebLink</i> sample tube
2050-M030-C1	TR30-C1	'Meinhard' Concentric glass nebulizer, Type C, <i>with PFA argon and sample connectors</i>
2050-M030-A3	WE02-4371	'Meinhard' Concentric glass nebulizer, Type A3, <i>with PFA argon and sample connectors</i>
2050-0710NT		High-Performance ' Tyfoon ' Concentric Glass Nebulizer, with <i>NebLink</i>
2050-4250		Boreal 'No-Break' HF-resistant Concentric Polymer Nebulizer
2550-2005		HF-resistant MicroFlow PFA Concentric Nebulizer, 360 ul/min uptake version
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0710		Burgener Teflon MiraMist Parallel-Path Nebulizer
3150-0720		Burgener PEEK MiraMist Parallel-Path Nebulizer
3150-5000	AM-5000	Burgener AriMist
3150-0730		Teflon T2100 'Work-horse' Parallel-Path Nebulizer
2050-0920		HF-resistant Fixed Crossflow Nebulizer, for Scott type chamber
2052-0300		PTFE Endcap for CG Nebulizer, 35/6 mm adaptor
2030-029L	N077-6052	'Althea' Cyclonic Spray Chamber, LHE, Optima 2000/4000, volume 50 ml.
2030-030L	N077-6053	'Libra' Cyclonic Spray Chamber Baffled-LHE, Optima 2000/4000, vol. 50 ml.
2030-0302		'Libra' , Baffled cyclonic glass Spray Chamber, with O-ring/glass nebulizer port
2030-0303		Jacketed 'Libra' , Baffled cyclonic glass Spray Chamber, vol. 50 ml.
2030-0307		'Electron' Mini-Cyclonic glass Spray Chamber, volume 20 ml.
2030-0309		'Positron' Baffled Mini-Cyclonic glass Spray Chamber, volume 20 ml.
2030-3157		HF-resistant 'ESI' Polypropylene Cyclonic Spray Chamber
2030-2929		HF-resistant 'Erebus' Cyclonic Spray Chamber
3304-5791C		Jacketed 'Althea' Cyclonic Glass Spray Chamber-LHE, for volatile samples
2030-0302L		Quartz 'Libra' Baffled Cyclonic Spray Chamber, Optima 2000/4000
3157-1101		Quartz 'Libra' Cyclonic Spray Chamber, O-ring free, fits in PC3 Peltier enclosure
EP-4157-1000		PC3 Peltier-cooled Cyclonic Spray Chamber, with 12-Volt Power Supply
EP-9100-5700		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2031-0160	N077-0338	Optima 4300DV / 7300DV Standard Quartz Torch, 1-slot, ES-1402-4000
2031-0170	N077-1500	Optima 4300DV / 7300DV Torch Body, 3-slots, Radial
2031-0180	N077-0344	Optima 2000/4000 Torch Body, Short version
2031-0190	N077-0343	Optima 2000/4000 Torch Body, No Slot
2031-0250	N077-5289	Optima 2000/4000 Torch Bonnet-DV, ES-1402-9001

2031-0400	N077-5227	Alumina Injector 0.8mm, Optima 4300
2031-0410	N077-5228	Alumina Injector 1.2mm, Optima 4300
2031-0420	N077-5177	Alumina Injector 2.0mm, Optima 4300
2031-0430	N077-0437	O-ring Kit for Optima 4300 torch
2031-0440	N077-0438	O-ring Kit for Optima 4300 Injector Assembly
2032-0700	0057-2615	Small Impact Bead 1/8"
2032-0540	N077-1531	Alumina Injector ID=1,2 Optima 4300V
2032-0740	N077-5014	Quartz Injector ID=2.0, L=141mm, Optima 2000/4000
2032-0750	N077-5224	Quartz Injector ID=3.0, L=141mm, Optima 2000/4000
2032-0760	N077-5223	Quartz Injector ID=1.6, L=141mm, Optima 2000/4000
2032-0770	N077-5226	Quartz Injector ID=1.2, L=141mm, Optima 2000/4000
2032-0780	N077-5225	Quartz Injector ID=0.8, L=141mm, Optima 2000/4000
2032-0800	N077-0322	Radial Purge Window, Optima 4000
201K-0700	N077-1116	Axial Purge Window, Optima 4000
2032-1503	ES-1503-8181	HF-resistant Semi-Demountable Torch with Sapphire Injector, ID = 1,8 mm.
2032-1513	ES-1513-8201	HF-resistant Semi-Demountable Torch with Platinum Injector, ID = 2,0 mm.

EPOND P/N°		Sample Introduction Accessory description
SPIRO-TMD	EP-4599-3002	' Spiro-TMD ' autonomous Membrane Separator, for solvent removal, with Power Supply
APEX-Q	EP-4857	APEX-Q high-sensitivity desolvation system
2019-4260		' XENA ' micro-flow meter / aspiration monitor, range : 5 - 1999 ul/min., mains powered
QuDIN-1011PE		Mermet-Todoli QDIN adaptor, with Peri-Pump tubes and connectors
2045-0700		' SOLVIT ' Argon Humidifier, with tubing and connectors
8100-2051		' PEARL ' Cleaning Device for concentric glass nebulizers - <i>suits all types</i>
6808-2280		' ISOBAR ' Internal Standard Kit, HF-resistant, with tubing and connectors

EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
EPS-00xx		2-tag Silicone P-P tubes, pkt 6, all internal diameters available
See web site		<i>Cone-terminated</i> Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40
DLK-075		<i>DrainLink</i> nebulizer adaptor and tube with 0,75mm ID, pkt 5
NLK-030		<i>Neb-Link</i> sample tubes, packet of 10, ID=0,30 mm, L=70cm.

Autosamplers and ICP / AAS Standard and Reference Solutions		
SC-8000-XXXX		Complete range of ESI SC-series autosamplers and parts - send for prices and service
ICP-standards		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards
		<i>... custom standards supplied - according to your specifications ...</i>



www.e-pond.biz