

2012		EPOND Accessories List - Thermo-Fisher ICP-OES instruments
Epond N°	OEM P/No.	ICAP-6300, 6500, IRIS Intrepid, Advantage, Axial + DUO
2050-0800N	13436500	Concentric Glass Nebulizer, K-type, <i>with NebLink sample tube</i>
2050-0800Q		Concentric Quartz Nebulizer, K-type, <i>with NebLink sample tube</i>
2050-0710N	8423 120 51431	Standard Concentric Glass Nebulizer, for iCAP series
2050-0710NT	14186300	Tyfoon High-performance concentric glass nebulizer, with NebLink capillary
2050-0650QD	8423 120 51331	Aerosalt Concentric Glass Nebulizer for iCAP series
2050-M030-K2	13960000	'Meinhard' Concentric glass nebulizer, Type K, 2 ml/min.
2050-4250		'NB-Class' (No-Break) Boreal Concentric Polymer Nebulizer, medium uptake version
2050-0920	12698900	TJA Fixed Cross Flow Nebulizer
2051-0100QD	8423 120 51321	V-Groove nebulizer, for iCAP series
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2042	ES-2042	Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041	ES-2041	Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
3150-0730	13638800	Burgener PTFE T2002 'High-Solids' Parallel-Path Nebulizer
3150-0710	8423 120 51471	Teflon MiraMist Burgener Parallel-Path Nebulizer, for iCAP series
BRC-203X		Polyethylene capillaries for Burgener nebulizer, pkt 10
BRG-204		Upchurch gas adaptor and argon tube for Burgener nebulizer
BRC-20		Upchurch gas adaptor for Burgener nebulizer
2052-0300		Teflon Spray Chamber End Cap, for Open-End Spray Chambers, 6mm nebs
2027-5150	13462000	TJA Axial / DUO Trace Cyclonic Glass Spray Chamber, 12mm drain
2027-515NB		TJA Axial / DUO, Glass Spray Chamber like 2027-5150, No Baffle
2027-5160	13612000	TJA Axial / DUO 'Althea' Cyclonic Glass Spray Chamber, with elbow
2027-5163		TJA Axial / DUO 'Libra' Baffled Cyclonic Spray Chamber
2027-5166		TJA Axial / DUO 'Electron' Mini-Cyclonic Glass Spray Chamber
2027-5169		TJA Axial / DUO 'Positron' Baffled Mini-Cyclonic Spray Chamber
2028-5150	8423 120 51411	Glass Cyclonic Spray Chamber, iCAP 6000 series
2028-5130	8423 120 51311	Glass Cyclonic Spray Chamber with baffle, iCAP 6000 series
2028-5135	8423 120 51441	Jacketed Glass Cyclonic Spray Chamber with baffle, iCAP 6000 series
2028-5260	8423 120 51251	Spray Chamber Adaptor for iCAP 6000 DUO series
2027-2929		EPOND 'Erebus' HF-resistant Cyclonic spray chamber
EP-4168-1000	ES-4168-1000	PC3 Peltier-cooled Cyclonic Quartz Spray Chamber, with nebulizer and Power Supply
EP-9100-6800		PC3-LT (low-temp.) Peltier-cooled Quartz Spray Chamber, reaches -20 degrees C.
2025-0230	13768900	TJA DUO Quartz Sleeve with Optical Port, ID = 22,5 mm minimum
2027-5100	126432-03	TJA Axial, Trace Torch with Quartz Aqueous Injector, ID = 1,5mm
2027-510SI	126440-03	TJA Axial, Trace Torch without Injector, (EP-1504-0200)
2027-510DV	138562-00	TJA DUO, Trace DV Torch with 1.5 Aqueous Injector
2027-510DVSI	126440-05	TJA DUO, Trace DV Torch without Injector
2027-510HF	138566-03	TJA Axial, Trace Torch 1.5 Alumina Injector
EP-1423-0100		Quartz 'Drop-In' Injector with ID = 1,0 mm
EP-1423-0150		Quartz 'Drop-In' Injector with ID = 1,5 mm
EP-1403-0100		Sapphire 'Drop-In' Injector with ID = 1,0 mm

EP-1403-0180		Sapphire 'Drop-In' Injector with ID = 1,8 mm
EP-1413-0180		Platinum 'Drop-In' Injector with ID = 1,8 mm
EP-1413-0200		Platinum 'Drop-In' Injector with ID = 2,0 mm
2025-0400	125407-01	TJA Radial, Axial and DV Quartz 1.0 Injector for Organics
2025-0410	125407-00	TJA Radial, Axial and DV Quartz 1.5 Injector for Aqueous
2025-0420	136502-01	TJA Radial, Axial and DV Quartz 2.0 Injector
2025-0430	135602-02	TJA Radial, Axial and DV Quartz 2.4 Injector
2027-5250	085207-00	TJA Torch Anti Arc Cap
2027-5260	085207-01	TJA Torch Anti Arc Cap / Flat Edge
2028-5100	8423 120 51241	iCAP 6000 series DUO Quartz Torch, older type
2028-5400	8423 120 51341	iCAP 6000 series Centre Tube, 1.0mm
2028-5410	8423 120 51351	iCAP 6000 series Centre Tube, 1.5mm
2028-5420	8423 120 51371	iCAP 6000 series Centre Tube, 2.0mm, L = 97mm
2028-5270	8423 120 51591	Duo Periscope Quartz Window/Sleeve
2028-5120	8423 120 51841	iCAP 6000 series DUO 'EMT' Quartz Torch, ML128050
2028-8072	8423 120 51971	iCAP 6000 series Centre Tube for EMT torch, 2.0mm, ML128072, L = 104,2mm
EPOND No.		Accessory Description
QuDIN-1011		Mermet-Todoli QDIN adaptor, with Peri-Pump tubes and connectors
APEX		APEX Desolvator, with PolyPro-ST nebulizer and SEDNA micro - flow meter
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply
2019-4260		XENA volumetric micro-flow meter / aspiration monitor, range : 10 – 1999 ul/min.
8100-2051		PEARL cleaning device - suits ALL types of concentric glass nebulizers
6808-2280		'ISOBAR' Internal Standard Kit, HF-resistant, with tubing and connectors
2045-0700		'SOLVIT' Argon Humidifier, with tubing and connectors
EPME-PG41-100		100-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm
EPME-PG41-500		500-ml bottle multi-element standard solution for ICP analyses, 28 elements @ 100ppm
		reference solutions made to measure - on request
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
ERT-040		EPOND 'RubeTubes' , Tygon-core sample tube connectors, packet of 40
PRT-040		EPOND 'RubeTubes' , PharMed sample tube connectors for organics, packet of 40
DLK-080		<i>DrainLink</i> nebulizer adaptor and tube, pkt 5
NLK-080		<i>Neb-Link</i> sample tubes, packet of 10, L=60cm.
		Autosamplers and ICP / AAS Standard and Reference Solutions
SC-8000-XXXX		Complete range of ESI SC-series autosamplers - please ask for details
ICP-standards		Certified ICP / AAS Standard solutions - http://www.e-pond.biz/standards
		reference solutions made to measure - on request
		 
		www.e-pond.biz