

2012		VG PQ series, Axiom, ExCell, PlasmaTrace
EPOUND N°	OEM P/No.	Part Descriptions
2050-0710N	1040894	EPOUND ' Caliber ' standard Concentric Glass Nebulizer, with NebLink sample tube
2050-0770Q		EPOUND A-type Concentric Quartz Nebulizer, with NebLink sample tube
2050-M030-C1	TR30-C1	' Meinhard ' Concentric glass nebulizer, Type C
2050-4151		' Zefyr ' 'NB-Class' (No-Break) Concentric Polymer Nebulizer, low-uptake version
2050-4251		' Boreal ' 'NB-Class' (No-Break) Concentric Polymer Nebulizer, medium uptake version
2550-2020		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 20 ul/min uptake version
2550-2000		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 50 ul/min uptake version
2550-2002		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 100 ul/min uptake version
2550-2003		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 200 ul/min uptake version
2550-2005		HF-resistant <i>MicroFlow</i> PFA Concentric Nebulizer, 400 ul/min uptake version
2550-2028		PolyPro-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2550-2024		PFA-ST HF-resistant concentric Nebulizer, with detachable sample uptake tube
2030-2028	ES-2030-28	PFA-HS HF-res. High-Solids concentric neb, with detachable sample uptake tube
2550-2045		Detachable sample capillary tube, ID = 0,15mm., nominal uptake 20 ul/min.
2550-2043		Detachable sample capillary tube, ID = 0,20mm., nominal uptake 50 ul/min.
2550-2042		Detachable sample capillary tube, ID = 0,25mm., uptake 100 ul/min.
2550-2041		Detachable sample capillary tube, ID = 0,50mm., uptake 400 ul/min.
		Capillaries with integrated probes - available for all common autosamplers
2052-0300		PTFE Endcap for CG Nebulizer, 35/6 mm adaptor
2045-1136		Standard ' Electron ' 20ml. Mini-cyclonic glass spray chamber
2045-1137		Standard ' Positron ' baffled Mini-cyclonic glass spray chamber
2045-1136Q		Quartz version of ' Electron ' Mini-cyclonic spray chamber
2045-1132		Standard ' Althea ' 50ml. Cyclonic glass spray chamber
2045-1132Q		Quartz version of ' Althea ' 50ml. Cyclonic spray chamber
2045-1142		Jacketed ' Althea ' 50ml. Cyclonic glass spray chamber
2045-1134		' Libra ' Baffled Cyclonic glass spray chamber
2045-1144		Jacketed ' Libra ' Baffled Cyclonic glass spray chamber
2045-0220	3204313	VG Jacketed Spray Chamber with Impact Sphere
2045-022Q		VG Quartz Jacketed Spray Chamber, with Impact Sphere
2045-0130	3600170	VG Quartz Conical Spray Chamber, for Peltier system
2045-0140		VG Horizon Conical Spray Chamber for 6-mm OD Nebulizers
2045-0160	3200841	Quartz Coaxial VG Spray Chamber, Water Cooled
2045-0170	3203739	VG Thermo Controlled Maxim Spray Chamber
2045-2929		HF-resistant ' Erebus ' Cyclonic Spray Chamber
2045-0800		HF-resistant Teflon Baffled Spray Chamber
EP-2124-0351	ES-2124-0351	35mm PFA spray chamber and PureCap End Cap for PQ2 and PQ3
EP-6124-0001	ES-6124-0001	Spray Chamber Cooling Block, anodized Al, plus Mounting Bracket No. 6124-1002
EP-6124-1002	ES-6124-1002	Mounting Bracket for attachment of Al cooling block to PQ2
EP-6124-2002	ES-6124-2002	Mounting Bracket for attachment of Al cooling block to VG Axiom
2055-0180	1040950	JC-13 Metal Spring Clamp
2045-0050	3201192	VG Standard Quartz Torch 1.5mm Injector
2045-0060	3201732	VG Standard Quartz Torch 2mm Injector
2045-0080	1200080	VG Demountable Torch Body, for 2045-0100
2045-0100		VG Semi-Demountable Torch, with Qtz Injector
2045-0090		VG Quartz Injector, for 2045-0100
6608-0252		VG Semi-Demountable Torch with Inlet elbow and Quartz Injector
6608-0244		Quartz Body for SDT 6608-0252
1200-0400		Semi-Demountable Torch with O-ring free Quartz Injector, ID = 1,5 mm
1200-0600		HF-resistant Semi-Demountable Torch with O-ring free sapphire Injector, ID = 1,8 mm
1200-0700		HF-resistant Semi-Demountable Torch with O-ring free platinum Injector, ID = 1,8 mm
1233-0150		Quartz Injector with O-ring free PFA base, ID = 1,5 mm
1203-0180		Sapphire Injector with O-ring free PFA base, ID = 1,8 mm
1213-0180		Platinum Injector with O-ring free PFA base, ID = 1,8 mm
1202-0200		Outer Quartz Torch for O-ring free Injectors
1202-9000		VG Axiom Shield Torch Bonnet
2045-0120	3203095	VG Torch Body CR Horizon
2045-0250		Standard VG Torch Bonnet, for PQ series, Micromass
2045-0260	1200274	VG Torch Bonnet Cap, Outside Diameter = 52.4mm
2045-0270	3203025	VG RF Insulator Shield, Outside Diameter = 140mm
2045-0280	3204703	VG Quartz Bonnet Axiom, O.D. = 72mm

2045-0290	3203789	VG Torch Shield, for PlasmaScreen option, suits VG PQ ExCell
2045-0300	3600169	VG Quartz Elbow
2045-0310		VG Elbow, Ball/Socket Pyrex
2045-0320	3200727	VG Elbow, Ball/Socket Quartz
2045-0360	3007383	VG UV Laser Sample Cell
2045-0400	3201104	VG Probe UV Laser Sample Cell
2045-0420		VG Laser Sample Cell, Micro Probe
2045-0430	3202287	VG Quartz Sleeve
2045-0440	3201759	VG Quartz Tube
2045-0450	3200176	VG Quartz Sleeve #2
2045-0460	3201070	VG Quartz Sleeve, Silica
2045-0470	3201568	VG Torch Tube Adapter
2045-0490	3200004	VG Quartz Socket Tube Adapter
EPOND No.	OEM P/No.	Cones – PQ series, PlasmaTrace, ExCell
VG1001-Ni	3004661	Nickel Sampler (P/No. 4001690 for pack of 3 **)
VG1002-Ni	3004866	Nickel Mini-Skimmer Cone
VG1004-Ni	3200860	Nickel Micro-Skimmer Cone (enhanced interface)
		<i>Cones above are available in aluminium, at nickel cone prices</i>
VG1006-Pt	3200262	Platinum Sampler Cone (15,2 mm Pt insert)
VG1006A-Pt	-	Platinum Sampler Cone (18 mm Pt insert)
VG1007-Pt	-	Platinum Mini-Skimmer Cone
VG1008-Pt	3201101	Platinum Micro skimmer Cone (Enhanced Interface)
VG1009-Pt	3202005	Platinum PT1 Skimmer Cone
VG1010-Ni	3200456	Nickel PT 1 Skimmer Cone
VG1011-Pt	3201941	Platinum PT2 Skimmer Cone
VG1012-Ni	3201585	Nickel PQE Skimmer Cone
VG1014-Ni	3201457	Nickel PT2 Micro Skimmer Cone
VG5001	3004382	Graphite Gasket for Sampler cones
VG5001A	3004382	Graphite Sampler Cone Gasket – old interface with screw holes
EPOND No.	OEM P/No.	Xi – series cones
VG1021-Ni-e	-	eSPEC Nickel Sampler Cone, 100% Nickel
VG1021-Ni/Cu-e	3600812	eSPEC Nickel Sampler Cone, Cu core
VG1022-Ni-e	3600811	eSPEC Nickel Skimmer Cone
VG1023-Cu	3600813	Copper Skimmer Adaptor
VG1026-Pt	1600335	Platinum Sampler Cone
VG-1028-Pt	N.A.	Platinum Skimmer Cone
2085-2206-05		Diamond paste for polishing and cleaning cones, 3-micron, 5-ml syringe
2085-2206-18		Diamond paste for polishing and cleaning cones, 3-micron, 18-ml syringe
EPOND No.	OEM P/No.	Description of ICP-MS Detector
EPM14562		Detector for PQ1 and PQ2 (with 12-12 rods)
EPM14562A		Detector for PQ1 and PQ2 (with SXP rods)
EPM4870V	2027814	Burle Channeltron detector – replacement now by Type 401A
EPM14562A		Detector for PQ3, ExCell (sequential)
EPM14214	AF 214	Detector for PQ3, ExCell (simultaneous)
EPOND No.		Accessory Description
APEX-Q-VG		APEX-Q Desolvator, with 2 PFA concentric nebs and 'Sedna' flow monitor
EP-4599-3002	ES-4599-300A	Spiro TMD Membrane Desolvation Module, with Power Supply
2019-4260		'XENA' micro-flow meter / aspiration monitor, range : 5 – 1999 ul/min., AC powered
8100-2051		'PEARL' Cleaning Device for concentric glass nebulizers
6808-2280		'ISOBAR' Internal Standard Kit, HF-resistant, with tubing and connectors
EPWP-0710	WP-710	Water Purification System - produces up to 30 litres/hour of UP water
		reference solutions made to measure - on request
EPOND N°	OEM P/No.	Description, tubes
See web site		All tube Internal Diameters available, 2 tag Tygon P-P tubes, pkt 12
See web site		Cone-terminated Peri-Pump tubes, for low-uptake apps., pkt of 10
ERT-040		EPOND 'RubeTubes', sample tube connectors, packet of 40
DLK-075		DrainLink nebulizer adaptor and tube with 0,75mm ID, pkt 5
ICP-standards		Certified ICP / AAS Standard solutions - www.e-pond.biz/standards
		reference solutions made to measure - on request
		www.e-pond.biz