

# PFA PureTip Pipette Tips

Often overlooked as a source of sample contamination the pipette tip can come into direct contact with the sample or acids that are added to the sample.

Only PureTips can deliver contamination free concentrated chemicals to the sample. PureTips can be reused hundreds of times and will not leach contaminants throughout their lifetime.

- Molded from high purity PFA fluoropolymer, minimizing contamination
- Chemically resistant to virtually all sample types, ideal for strong acids, alkalis, organics and biologicals
- Long life, re-usable
- Use with standard Eppendorf type 1000 $\mu$ L & 100 $\mu$ L pipetters



PureTip Pipette Tips are available in 1000 $\mu$ L & 100 $\mu$ L sizes and can be used with standard Eppendorf type pipetters



ES-7000-1011 - 1000  $\mu$ l Pipette Tips  
Quantity of 60, pre-cleaned



ES-7000-0111 - 100  $\mu$ l Pipette Tips  
Quantity of 96, pre-cleaned