

SC-14 CLEAN ENCLOSURE

The Clean Enclosure has been developed for the SC-14 high capacity autosampler. This enclosure helps to maintain sample integrity while samples await analysis.

- ✓ Protects samples from airborne contamination.
- ✓ Clear acrylic construction, provides sturdy protection to the SC-14.
- ✓ Covers only the sample platform
- ✓ 2 large door openings for easy access to the sample racks.
- ✓ Clean Enclosure doors easily open and lock away closed for rapid access to samples.

P/N SC-1107
SC-14 Enclosure



SC-14 autosampler with clean enclosure attached



Enclosure doors easily open for access to sample racks