

SC-2 CLEAN ENCLOSURE

The Clean Enclosure has been developed for the SC-2 ultraclean autosampler. This enclosure helps to maintain sample integrity while samples await analysis.

- ✓ Protects samples from airborne contamination.
- ✓ Clear acrylic construction, provides sturdy protection to the SC-2.
- ✓ Covers only the sample platform
- ✓ Large door opening for easy access to the sample racks.
- ✓ Clean Enclosure door easily opens and locks away closed for rapid access to samples.

P/N ES-8199-1207
SC-2 Enclosure



SC-2 autosampler with clean enclosure attached



Enclosure doors easily open for access to sample racks