


NEW MP² micro peristaltic pump

- Very Low Pulsation
 - Precision/internal standard addition at low-flow rates.
 - Faster analysis without losing precision.
- Wide Dynamic Range
($< 1\mu\text{L}/\text{min}$ - $> 10\text{mL}/\text{min}$)
- Optimized for Low Flow Rates
- Stand Alone and Integrated Versions
- One to Six Channels
- Easy to Use
 - Aligned Pressure Shoes
 - Easy Tube Replacements
 - Controlled by SC-Software




APPLICATIONS:


- Precise On-line Internal Standard Addition




- Low-Flow Applications




- Microchromatography



- High Precision - SSI



ESI Pump - $30\mu\text{L}/\text{min}$, orange/green tubing (0.38mm i.d). 1.4% RSD



Standard Pump - $30\mu\text{L}/\text{min}$, orange/green tubing (0.38mm i.d). 3.9% RSD